Listen to your Heart:

Heartbeat Sound Segmentation & Classification

September 11, 2019

Boikanyo Radiokana & Elias Sepuru

School of Electrical & Information Engineering
University of the Witwatersrand
South Africa



Agenda



Introduction

Objectives


Background
Heartbeat Sounds Categories
Related Work
Project Setting

Methodology Data Acquisition



 CVDs are the leading causes of death globally - WHO.

Annual Number of Deaths by Cause





- CVDs are the leading causes of death globally - WHO.
- Currently used method to check for CVDs is Cardiac Auscultation (CA).





- CVDs are the leading causes of death globally - WHO.
- Currently used method to check for CVDs is Cardiac Auscultation (CA).
- ► CA is a difficult skill to acquire.



Correct diagnosis using CA in USA, Canada & UK respectively.




CVDs are the leading causes of death globally - WHO.

- Currently used method to check for CVDs is Cardiac Auscultation (CA).
- ► CA is a difficult skill to acquire.
- ► People are not aware of their heart conditions.

Awarness of Heart Condition

America - 2016



- Know Their Heart Condition
- Don't Know Their Heart Condition



Easily accessible & reliable heart diagnosis systems would help reduce deaths due to CVDs.

Objectives

➤ To segment Heartbeat sounds (HSs) based on the location of S1 (lub) S2 (dub) in Normal HSs.




Objectives |

➤ To segment Heartbeat sounds (HSs) based on the location of S1 (lub) S2 (dub) in Normal HSs.



► Create models that will enable preliminary screening of CVDs





This project deals with classifying HSs into the following categories:

- 1. Normal HSs
- 2. Murmur HSs
- 3. Extra Heartsounds
- 4. Exrasystole HSs
- 5. Artifact



Normal HSs

lub...dub.....lub...dub....





Murmur HSs

```
lub...***..dub.....lub...***..dub......

or

lub....dub...***...lub....dub...***...
```





Extra HS

```
lub.lub...dub.....lub.lub...dub.....


or

lub...dub.dub.....lub...dub.dub.....
```






Extrasytole HSs



Background Heartbeat Sounds Categories




Artifact Sound Not an actual HSs.



Background Heartbeat Sounds Categories



Can you guess the categories?





Strunic's attempt to classify HSs with ANN.

85±7.4%

Accuracy when classifying simulated HSs with no noise.



Accuracy when classifying real life HSs with noise.

Background Project Setting




To make this project applicable to real world situations, two datasets recorded in real life settings will be used. Both datasets contain excessive background noise.

Methodology Data Acquisition



Dataset A




- Recorded by the general public
- Device iStethoscope Pro lphone app
- ► Sampling Freq 44100Hz
- Contains excessive background noise

Methodology Data Acquisition



Dataset B



- Recorded from a hospital by Medical Practitioners
- ► Device Digital Stethoscope
- ► Sampling Freq 4000Hz
- Contains background noise